
TECNIKA

1- 3	Mode d'emploi	F	S Användarmanual	14-16
3- 5	User Manual	GB	FIN Käyttäjän opaskirja	17-19
5- 7	Bedienungsanleitung	D	DK Brugermanual	19-21
8-10	Modo di uso	I	NL Gebruiksaanwijzing	21-24
10-12	Manual de utilização	P	GR ΟΔΗΓΙΕΣ ΣΡΗΣΗΣ	24-26
12-14	Manual para el usuario	E		

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TECNICA

1. NOTICES

⚠ This symbol, placed on the rear label warns the user to refer to the manual.

TECNICA is an endodontic micromotor suitable with every kind of Ni-Ti instruments. Its usage is allowed only in hospital environment, in clinics and in dental offices. It must be used only by specialized personnel in dentistry. **The manufacturer accepts no responsibility in the event of:** a) actions or repair carried out by person not authorized by the manufacturer or importer. b) a use on an electrical system not in conformance with the terms of regulation CEI 64-4. c) a different use from the one specified on this booklet.

2. COMPONENTS

1) CENTRAL UNIT Ref. : TCK/JUNT 2) PEDAL Ref. : TCK/PED 3) MICROMOTOR 30V 20W: TCK/MIC 4) CONTRANGLE SUPPORT Ref. :TCK/SPP 5) PROTECTION CAP FOR MICROMOTOR Mod ATR/PRT 6) " ATR SOFTWARE UPDATE: Tecnica ATR software 2000

3. DESCRIPTION

FRONT PANEL FIG. 1

1. "Motor" switches micromotor on and off.

2. "F+R" activates motor running CW/CCW. (Torque value is only 100). When is activated this function (LED = On), "A\Rev" is disarmed.

3. "A\Rev" Allows Automatic Reverse function turn On or Off (Available only if the set torque values are different from 100).

4. "Rev" Allows micromotor rotate CCW.

5. "system" Allows to change up to three different systems: LC = Large Canal ProFile Sequence, RC = Regular Canal ProFile Sequence (ProFile Sequence) GT = GT files Sequence.

6. "ctrl" Turns on or off some functions together with other keys, as follow:

• "ctrl"+"Rev" = turns on or off buzzer for CCW motor rotation.

• "ctrl"+"A\Rev" = turns on or off buzzer that, with reduced torque (different from 100), warns the user he is reaching the set value.

• "ctrl"+"reset" = sets programs 1,2,3 in all system at default values.

• "ctrl"+"reduct." = auto-check of contrangle. Please make this check monthly.

• "ctrl"+"F+R" = Enters menu to set time (sec/100) for CW (use torque. + and - key) and CCW (use rpm + and - key)

• "ctrl"+"Motor" = Enters menu to set time (sec/100) to restart motor in CW rotation during autoreverse without removing the foot from the pedal. Use rpm + and - key to set delay time.

• "ctrl"+"Torque Max" = Enters menu to set percent of increase or decrease for the buttons " Tmax" and " Tmin". (Default value = 25%)

7. Display easy reading LCD alphanumeric display.

8. "reduct." Set reduction rate of used contrangle (1:1 - 15:1 - 16:1 - 18:1 - 20:1)

9. "reset" Together with " Control " key sets programs 1,2 e 3 in all systems at default values.

10. "torque" Increases or decreases torque values from 0 (minimum) to 100 (Maximum).

11. "rpm" Set micromotor rotation speed. Max speed with reduction rate 1:1 is 12800 rpm.

12. "Torque Min" Decreases the torque value in use according with percent set in the menu "ctrl + Torque Max".

13. "Torque Max" Increases the torque value in use according with percent set in the menu "ctrl + Torque Max".

14. "3" Program n° 3. (See after)


15. "2" Program n° 2. (See after)

16. "1" Program n° 1. (See after)

FRONT PANEL Fig2

1. Pedal connector
2. Micromotor connector

REAR PANEL Fig 3

1. Main switch
2. Change voltage (110/220 V).
3. Fuses.
4. Main supply
5. Serial input (J1). *

*Use only with " ATR SOFTWARE UPDATE: Tecnica ATR software 2000 "

This serial port (j1) allows to connect TECNICA to a PC in order to update the software. Directions for use: see " ATR SOFTWARE UPDATE: Tecnica ATR software 2000 "

4. INSTALLATION

The unit in its original packaging can be stored at temperature from 5° C (41°F) to 65° C (149°F) with relative humidity lower than 90%.

Open the box cutting adhesive tape. Remove the unit and its accessories in their original packaging and put on a flat surface. Remove unit and accessories from packaging cutting their protection. Let at least 150 mm free space around the unit for ventilation.

Before installing make sure that unit has no damages. Possible damages caused by carrying must be reported to the dealer within 24 hours from receiving. Verify the content of the box according to paragraph 2. COMPONENTS.

** This unit must be used only by specialized person. Do not allow no authorized people to touch the unit or play with it and to use it for any other purpose.

** Unit must not be placed near direct or indirect heat fonts.

** Unit must not be placed in a humid environment or in contact with any kind of liquids.

- ** Let air circulate around the Unit for its cooling.
- ** It's forbidden and dangerous to modify Unit characteristics; in this case A.T.R. S.r.l. will not be responsible for possible damages.

Environmental conditions:

1. Altitude: lower than 3000 m (10000 ft)
2. Temperature: 18°C. / 40°C (64° F / 140° F)
3. Relative Humidity: < 80%
4. Current fluctuations cannot exceed +/- 10% nominal voltage
5. On the motor can be attached only CE conrtriangle or legally sold conrtriangle.
6. Do not introduce any other item into the micromotor.

Not suitable for use in the presence of flammable anesthetics or oxygen.
Micromotor duty cycle: 30 seconds on / 10 seconds off.

5. USE

Before switching the unit on, be sure that voltage on rear label is the same of the main supply. Before connecting unit to main supply complete mounting of all its parts. Earthing the unit is obligatory by law. Connect **TECNIKA** to a valid earthing. No manufacturer responsibility for damages to things or people due to non observance of this direction. **Connector J1 on the rear panel must not be connected with any device. Directions to use this port will be given with updating kit provided by ATR.** The micromotor must work in intermittent use (interval of minimum 10 seconds every 30 seconds of use) to avoid micromotor overheating. When completed all the connections :

a) Push main switch. (Fig. 3.1)

At first lighting, the Unit shows on display for 3 seconds **"Welcome in Digital Endodontics"**.

Display – Last fixed system: RC1* T:7 16:1 250
 PF04 25/20 Tmin:/15

First row (from left to right) shows: • RC: System in use • 1: Program in use (1,2,3) • : "Special Autoreverse" in use (* on) • T7: Torque value • 16:1: Contra-angle reduction rate • 250: Speed.

Second row shows: PF04: *Type of file (Profile 04 06, GT Rotary File., Orifice Shapers) and code number of file to be used * If display shows "Tmax" or "Tmin" it means that to obtain the suggested torque for type and code of instruments showed on display after this symbol is necessary to push the button like the showed symbol. With this operation torque value will be modified with a pre-selected percent (Default = 25 plus or minus). In this situation the torque value showed on the first row will blink. To come back to the previous torque value will be sufficient to push the previous button again.

"SYSTEM" button selects: "RC" = Regular Canal Preparation; "LC" = Large Canal Preparation; "GT" = GT files • **"PROGRAMS"** 1, 2, 3 buttons select programs of each kind of instrument to be used with each system • **"REDUCTION RATES"** available 4 reduction rates: 15:1,16:1,18:1,20:1. (1:1 = reference scale.)

"TORQUE VALUE" : with reduction rates 15:1,16:1,18:1,20:1 torque available at instruments mounted on conrtriangle is always the same, because **TECNIKA** calculates the reduction rate set with the button "reduct.". With 1:1 rate torque is the same one of the motor shaft.

"SPEED" : with reduction rates 15:1,16:1,18:1,20:1 speed available at instruments mounted on conrtriangle is always the same (differences around 2% due to the calculus approximation), because **TECNIKA** calculates the reduction rate set with the button "reduct.". With 1:1 rate speed is from 1600 to 12800 rpm.

All programs (1, 2, 3, [see after]) have default values.

N.B. Default values of speed and torque must be considered "Suggested Values."

-Profile 04/06, GT Rotary Files, Orifice Shapers are Dentsply Maillefer registered Trade Marks.

" F+R" Function

Whenever activate the key "F+R", motor rotates with alternative movement forward/reverse . No torque control is allowed with F+R movement. (T = 100) Moreover it is possible to choose time of rotation forward and time of rotation in reverse from minimum of 0.05 sec to maximum of 2.5 sec. To do this:

- Push " CTRL" key and "F+R" key together
- By " TORQUE" key increase or decrease time for forward
- By " RPM" key increase or decrease time for reverse
- By " T Min " key save data and esc.

"Autoreverse" Function

Tecnika has two different functions for "Autoreverse": ("Arev " key = yellow LED on).

Function 1 or " Normal Autoreverse":

Whenever the value of set torque is reached, the motor automatically begins the rotation in reverse. To return in the forward rotation, it is necessary to lift the foot from the pedal: motor stops. When the pedal is pressed again the motor begin to turn forward.

Function 2 or " Special Autoreverse":

Whenever the value of set torque is reached, the motor automatically begins the rotation in reverse and automatically begins to turn forward after a time selected by the operator. To select " Normal Autoreverse" or "Special Autoreverse" do as follows:

- Push "CTRL" key together with "MOTOR" key
- By "RPM" key increase or decrease time for "Restart Delay"
- If "Restart Delay" time is 0.00, " Normal Autoreverse" activates.
- If "Restart Delay" time is different from 0.00 "Special Autoreverse" is activated. "Restart Delay" can be increased up to 2.5 sec.
- By " T Max" key save and "T Min" key esc.

NB1: When "Special Autoreverse" function is on an asterisk appears on display.

NB2: When "Autoreverse" is off ("Arev" key = yellow LED off) and set torque is reached, motor stops. It will be necessary now to lift foot from pedal, manually set "Reverse" function ("REV" key = yellow LED on) and press pedal again. To enter forward rotation: lift foot from pedal and manually disconnect "Reverse" function ("REV" key = Yellow led off).

" Torque Max" and " Torque Min" Functions


"Torque Min" decreases the torque value in use according with percent set in the menu "ctrl + Torque Max".

"Torque Max" increases the torque value in use according with percent set in the menu "ctrl + Torque Max". 25% (plus or minus) is the default value. To change this value:

- Push "CTRL" key together with "Torque Max" key
- By "RPM" key increase or decrease the percent value of "Torque Max" (minimum 10% - Maximum 40% -)
- By "Torque" key increase or decrease the percent value of "Torque Min" (minimum 10% - Maximum 40% -)
- "Torque Min" Key esc.
- "Torque Max" key save and esc.

6. DEFAULT VALUES

In every programs of every systems reduction rate is:16:1 , rpm is 250 and "A/Rev" key is activated (LED = ON)

7. CLEANING

To clean and disinfect the external surfaces use a wet cloth.
The box is not watertight.

Don't dip any component in ultrasonic cleaners.

TECNIKA parts that can be introduced in autoclave (saturated steam 134°) • Micromotor • Micromotor's cable • Micromotor's connector • Micromotor support
Micromotor can not be used until its temperature goes under 40 °C.

It is obligatory to use the protection cap on the motor going into autoclave.

8. MAINTENANCE

Check periodically the integrity of main supply cord and if there is any damage replace with another with the same characteristics (220V 110V 500 mA T/ 1A T) approved by Notified Body or ask us for replacement.

FUSE REPLACEMENT

FUSE: 2 x 500mA T @ 220V or 2 x 1A T @ 110 V

• Switch off the unit • Disconnect the unit from the main power supply • Remove fuseholder cover • Substitute fuses with the same type on the label • Put on the fuseholder cover • Pay attention that on fuseholder's window appears the right voltage because has also function of voltage change.

9. TECHNICAL DATA

MECHANISCH DIMENSION: 165x190x105 • MATERIAL: BAYER FR 90 V0 • WEIGHT: 1.5 kg • NOISE: < 25 dBA
ELECTRICAL SUPPLY VOLTAGE: 220Vac 50Hz / 110Vac 60 Hz ± 10% • MAX INTENSITY: 0,2 A - MAX POWER: 25 VA • FUSE: 2 x 500 mA T@ 220 V 2 x 1A T@ 110 V CLASS / TYPE: 1 / BF, applied parts: micromotor contra angle and Ni-Ti file • CENTRAL UNIT: IP31 - Ordinary • MICROMOTOR: IP65 • PEDAL: IP21 - Ordinary

"IN CASE OF..."

1- The unit doesn't switch on: a) Verify that main cable is correctly inserted in the main socket and that the voltage corresponds to that written in fuseholder and verify the integrity of the fuses.

2- The footswitch doesn't activate the motor: a) Verify that the footswitch cable is correctly inserted in its receptacle b) Verify the operation of the motor by "motor" key on the keyboard of the central unit. - Perform the operation "CONTROL+ REDUCT".

3- The unit (or the display) doesn't work correctly: a) implement the "GENERAL SYSTEM RESET":

- switch the unit on.

- push "CONTROL" key and at the same time the "SYSTEM" key

- push "CONTROL" key and at the same time the "REV" key

- push "CONTROL" key and at the same time "3" key

- push "CONTROL" key and in the same time "SYSTEM" key

- push "CONTROL" key and in the same time "1" key

- push "1", "2" and "3" keys.

N.B.: Go on with this procedure only if it is strictly necessary because all working parameters inserted by operator will be lost and default parameters will be loaded .

b) Reload TECNIKA Default Software by "ATR Software 2000" (see user manual).

4- The motor doesn't start: a) Verify that the motor plug is correctly inserted in its receptacle • b) Verify the right fitting of the conrtriangle and possible damages on it • c) remove the conrtriangle, set the maximum speed of the motor (see user manual) and retry.

5- The motor lose power: • a) verify possible damages of the conrtriangle.

AFTER SALE SERVICE: ATR Hotline ☎+39 05 73 22325 • [http://: www.tecnikaonline.com](http://www.tecnikaonline.com)

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