# endoPocket

Endodontic electronic control system

GB

- Operating instructions •

**GB**

*Dear Customer,*

*Thank you for purchasing **ATR EndoPocket** endodontic electronic digital control system.*

*This system is very easy to use and it is equipped with advanced functions that have been studied and realised by the research and development department at Advanced Technology Research. This manual is provided for safety reasons and is therefore an integral part of the equipment. We recommend always keeping it within reach.*

*Since the technical documentation supplied is part of the product, it always must be kept within easy reach.*

*Please read these Operating Instructions carefully and familiarize yourself with the **ATR EndoPocket**. If you encounter any unsolvable problems after thoroughly studying the Operating Instructions, contact your dental depot.*

*Please do not hesitate to contact us concerning any doubts or problems related to this manual.*

**Contents**


1. Precautions	5
Disposal of the equipment	5
Notices	5
Adverse reactions	7
Precautions Regarding Electrostatic Discharge (ESD)	5
2. Part list	9
3. Description	10
Control and Operating Elements	11
Foot Pedal Control	12
Warnings	13
Environmental Conditions for Use	13
4. Installation	14
5. Use	15
6. Use with EP-Fill	19
7. Cleaning, disinfecting and sterilising	21
8. Maintenance	22
Servicing	22
Supplementary maintenance	22
9. Specifications	23
10. Warranty Conditions	24
11. Troubleshooting	25
12. Declaration of Conformity	26
13. Accompanying documents	27
Electromagnetic emissions	27
Electromagnetic immunity	28
Recommended separation distance	29
14. Authorized Service Center	31

## 1. Precautions

- ✓ The original packaging can be stored in environments whose temperatures range between + 5°C to + 65°C and whose relative humidity is less than 80%.
- ✓ Do not use in the presence of oxygen and flammable mixtures.
- ✓ Avoid acute angles folding to damage micromotor, footpedal and power cables.

Do not sterilize **ATR EndoPocket** central unit in autoclaves or ultrasonic cleaners for any reason.

### Disposal of the equipment

It generally applies that any disposal of this product must comply with the relevant national regulations. Please observe the regulations applicable in your country. Within the European Economic Community, Council Directive 2002/96/EU (WEEE) requires environmentally sound recycling/disposal of electrical and electronic devices.

Your product is marked with the adjacent symbol. Disposal of your product with domestic refuse is not compatible with the objectives of environmentally sound recycling/disposal. The black bar underneath the “garbage can” symbol means that it was put into circulation after Aug. 13, 2005. (see EN 50419:2005) Please note that this product is subject to Council Directive 2002/96/EU (WEEE) and the applicable national law of your country and must be recycled or disposed of in an environmentally sound manner. Please contact your dealer if final disposal of your product is required.

### Notices

*These warnings are provided for safety reasons. Please read them carefully before installing or using the equipment for the first time.*

*It is very important to make sure that this manual is stored together with the unit so that it will also be available at any time in the future. In case the unit is sold or the right to use it is transferred to a different owner, make sure that the manual always remains located near the unit. This is necessary to ensure that the new owner can use the unit properly and has all required information.*

*This symbol appears on the rating plate to remind the user to observe the manual.*

## 1. Precautions

**ATR EndoPocket** is an electrical digital control system for endodontics. The system complies with Directive **EEC 93/42 Medical Device**.

- ✓ It can only be used in suitable places (clinics and dentist practices) and by physicians and specialists regularly certified to perform dentistry.
- ✓ **ATR EndoPocket** requires special precautions concerning EMC and must be installed and commissioned in compliance with the EMC information contained in this instruction manual.
- ✓ To prevent any risk of electromagnetic interference, do not use electromedical devices or those of any other nature in proximity of **ATR EndoPocket**. The equipment emits electromagnetic radiation within the limitations set by current legislation (**EN 60601-1-2**).
- ✓ The equipment can be used in environmental conditions whose temperatures range from +18 °C to +40 °C and whose relative humidity is less than 80%.
- ✓ **ATR EndoPocket** does not require checks or preventive maintenance.
- ✓ Follow the manufacturer's instructions to use the files.
- ✓ **ATR EndoPocket** is not supplied with endodontic files and reducer contra-angles.
- ✓ Any contra-angle can be used with **ATR EndoPocket** as long as it is equipped with an ISO coupling in compliance with ISO Standard 3964.
- ✓ **Torque accuracy and rotation speed are only guaranteed when the following contra-angles are used: Sirona 6:1, ATR 6:1 and ATR 16:1. Torque and speed tolerance are equal to 10% with respect to the value indicated on the monitor. Torque and speed are the effective file torque and speed only if the user has selected the reduction ratio coinciding with the ratio of the contra-angle effectively used.**
- ✓ The maximum sterilisation temperature allowed for the cable and micromotor is 135°C. Before sterilising the motor, always remove the contra-angle and put on the safety cover.
- ✓ The voltage supply must fall between: 100 - 240 V, 47 - 63 Hz.
- ✓ In the event of malfunctioning, suspend work and contact the service centre.
- ✓ **Never insert an ATR EndoPocket memory stick into a PC USB port. Never plug a PC external memory (removable disc) into ATR EndoPocket's USB port. Never connect ATR EndoPocket to a PC using a female-female USB cable.**

## 1. Precautions

GB

The manufacturer declines all responsibility in the event of:

- a) use other than the ones specified in this instruction manual;
- b) interventions or repairs performed by people not authorised by the manufacturer or, on its behalf, by the importer;
- c) use with an electricity network not in compliance with the provisions of Standard **IEC 64-8 section 710**;
- d) use of different and/or non-original components (with respect to those specified in the paragraph PART LIST).

### ***Adverse reactions***

To date, no contra-indications or undesired effects deriving from normal use of **ATR EndoPocket** in a clinical setting have been reported.

The manufacturer **Advanced Technology Research S.r.l., Via S. Donato 1, 51100 - PISTOIA, ITALY** declines all responsibility for events deriving from clinical use of the equipment in question, and in particular in the case of breakage of the file in the root canal.

### ***Precautions Regarding Electrostatic Discharge (ESD)***

**ESD** is the acronym for **Electro Static Discharge**.

#### **WARNING**

**Do not touch the connector pins or make connections with these connectors unless protective measures against ESD have first been taken.**

Protective measures against ESD include:

- Procedures targeting limitation of electrostatic charges (e.g.: air conditioning, humidifying ambient air, installation of non-conductive floor coverings, use of non-synthetic clothing).

## 1. Precautions

- Discharging the operator's body charge to ground connections or large-size metal objects.
- Grounding oneself through an arm band.

We recommend all persons using this device be informed on the meaning of this warning and that adequate training is provided on the physics of electrostatic charges, which may build up in a dentist's office, and on the potential damage to electronic components in the event of contact by an electro-statically charged USER.

**ESD** is the acronym for **E**lectro **S**tatic **D**ischarge.

Any discharge implies a previous charge.

Electrical charges are generated whenever two bodies move against one another; for example, when walking (soles of shoes on the floor) or driving (tires on the road surface). The charge is greater in conditions of low (rather than high) ambient humidity and with synthetic rather than natural materials (clothing, floor coverings).

The following empirical rule provides an idea of the voltages involved in electrostatic discharges.

An electrostatic discharge:

- will be felt if it is 3,000 Volts or greater
- will be heard if it is 5,000 Volts or greater (bang, hiss)
- will be visible if it is 10,000 Volts or greater (flash).

The currents involved in these discharges are in the order of 10 amperes. These currents are innocuous to man because they only last a few nanoseconds. **1 nanosecond = 1 billionth of one second = 1/1,000,000,000 of one second.**

In the case of voltage differentials greater than 30,000 Volts per centimetre, the charge is compensated (electrostatic discharge, flash, discharge).

It is easy to understand how the integrated circuits permitting the miniaturisation of the electronic device in question are sensitive to electrostatic discharges.

Voltages that are imperceptible to the user can cause a current discharge sufficient to burn the microprocessor. Damage of the single integrated circuits can cause incorrect operation or the breakdown of the entire device.

**To prevent this, it is necessary that the user is electro-statically discharged before using the device.**

## 2. Part list

**GB**

**ATR EndoPocket** includes the following components:

1. Central unit
2. Micromotor with sterilizable cable (length 2.7 m) and connectors
3. Micromotor protection cover
4. Foot pedal control
5. External power supply, model Cincon Electronics Co. Ltd, TR30RAM240
6. USB memory stick
7. Micromotor support
8. Instruction manual

### 3. Description


**ATR EndoPocket** is an endodontic unit with a sophisticated design that takes up relatively little space, thanks to its foldable display.


### 3. Description

GB

#### Control and Operating Elements

1. Main switch

### 3. Description

2. **Function:** selects the file systems in use. Allows reloading the unit's default parameters (by simultaneously pressing Function ▲ and Function ▼).
3. **Files:** selects the file being used. Each file is associated with the correct torque and speed values.
4. **Torque:** modifies the torque of the file being used.
5. **Rpm:** modifies the rotation speed of the file being used.
6. **Reverse:** (LED on) sets a counter-clockwise rotation direction for the file. By holding down the button, saves the file parameters on the memory stick provided.
7. **A/Rev:** (LED on) sets the automatic reverse for the micromotor. The file changes rotation direction as soon as the torque set is reached. By holding down the button, reloads the file parameters from the memory stick provided.
8. **Cal:** calibrates the contra-angle to ensure torque accuracy each time the contra-angle is replaced or lubricated.
9. **Anatomy:** reduces the torque with respect to the torque set, depending on the difficulty of the endodontic treatment (root canal morphology).
10. **Red:** allows choosing between several reduction ratios of the contra-angle being used.
11. **Motor:** starts the motor from the keyboard rather than by using the pedal. The motor starts if it is pressed once. Pressing it again stops the motor.

#### **Foot Pedal Control**

**ATR EndoPocket's** micromotor can also be started up using the pedal. Contrary to the "Motor" button on the keyboard, the micromotor operates as long as the switch remains pressed. Pushing down on the pedal once after the motor has been started using the keyboard stops the motor (the pedal overrides the button).

### 3. Description

GB

#### **Warnings**

- ✓ Before installing the unit, ensure it has not been damaged. The retailer must be informed of any damage that occurred during transportation within 24 hours of receiving the goods.
- ✓ The unit must not be placed in humid environments or those in constant contact with liquids.
- ✓ The unit must not be placed close to sources of direct or indirect heat.
- ✓ Do not use cellular phones in proximity of **ATR EndoPocket**.
- ✓ This equipment must only be operated by specialised personnel. Any use different from its intended use is prohibited.
- ✓ It is strictly prohibited to alter or modify the characteristics of the product; the company **Advanced Technology Research S.r.l., Via S. Donato 1, 51100 - PISTOIA, ITALY** shall not be held liable in the event alterations or modifications are made to the product.
- ✓ The unit emits electromagnetic radiation within the limitations set by current legislation.


#### **Environmental Conditions for Use**

- ✓ Use: indoor
- ✓ Altitude: below 3,000 m
- ✓ Environment temperature: 18°C - 40°C
- ✓ Relative humidity: < 80%
- ✓ Voltage supply fluctuations must not exceed +/- 10% of the rated voltage (230 V).

## 4. Installation

- ✓ Open the box by cutting the adhesive tape. Avoid using sharp objects, which might damage the equipment.
- ✓ Delicately remove the unit and accessories from their packaging and place them on a level surface.
- ✓ Ensure all components are present as specified in the paragraph “Part list”.
- ✓ Ensure the serial number located on the label (on the lower part of the box on the external packaging) corresponds to the serial number located on the fiscal and transportation documents accompanying **ATR EndoPocket**.
- ✓ Plug the external power supply movable connector in the power connector located on the back of the unit (A). Ensure the connection is solid.

- ✓ Plug the micromotor connector in the connector on the front right side of the unit (red, pin 9) and the pedal connector in the unit's front left connector (blue, pin 4).
- ✓ The connectors have a spring retaining system preventing accidental release. To disconnect the cables, hold the connector in the central part and pull it out.

## 5. Use

GB

- ✓ Connect the external power supply to the socket. Power voltage and frequency must be between 100-240 V and 47-63 Hz.
- ✓ Remove the sterilisation safety cover from the micromotor.
- ✓ Turn the unit on by pressing the main ON switch 1. Once on, a welcome screen indicating the software release version will appear:

ATR EndoPocket  
[www.atrdental.com](http://www.atrdental.com)


- ✓ The micromotor must be operated **intermittently** (30 seconds ON, 10 seconds OFF) to prevent it from overheating.
- ✓ Connect the contra-angle to the micromotor and keep the **Cal** button pressed until the micromotor starts turning to measure the inertia of the contra-angle. Warning: during calibration, the micromotor varies its speed from the minimum value (1600 rpm to the motor shaft) to the maximum value (12800 rpm to the motor shaft). For safety reasons, we recommend performing calibration without an endodontic file connected to the contra-angle.

=====  
Software version x.x  
DD/MM/YYYY  
=====

## 5. Use


During the calibration process, the screen appears as follows:

Press any button to exit the calibration process. The screen appears as follows:

- ✓ Once ON, the unit displays the following screen:

## 5. Use

GB

- ✓ **First line:** indicates the *file* system in use. The system can be selected using the **Function** button.
- Second line:** indicates the file selected. The file can be selected using the **Files** button.
- Third line:** indicates the reduction ratio of the contra-angle (default: **6:1**) to the left and the file rotation speed (associated with the contra-angle selected) to the right. Reduction ratio and speed can be modified using the **Red** and **Rpm** buttons respectively.
- Fourth line:** indicates the stop/autoreverse torque for the file selected to the left and the morphological difficulty (**Anatomy**) to the right. Torque and Anatomy can be modified using the **Torque** and **Anatomy** buttons respectively.

**An:** the following are the values that can be selected for **Anatomy**:

- **E (default): easy.** The file stop/autoreverse torque is equal to the value set.
  - **M: medium.** The file stop/autoreverse torque is equal to the value set less 15%.
  - **D: difficult.** The file stop/autoreverse torque is equal to the value set less 30%.
- ✓ **In the event the user changes the reduction ratio, torque and speed on every file, the unit saves the changes. To return to the default parameters, simultaneously press Function ▲ and Function ▼. **WARNING: the changes are memorised if they remain constant for at least 3 seconds.****
- ✓ **A/Rev:** when the LED is on, the file inverts its rotation direction (from counter-clockwise to clockwise) as soon as the torque set on the display is reached. Lift your foot from the pedal and press down on the pedal again to return to counter-clockwise rotation. When the LED is off, the file stop as soon as the torque set is reached. **When the file reaches 75% of the torque set, the unit emits an acoustic signal.**  
Using the **A/Rev** button, torque, speed and reduction ratio parameters can be loaded from the memory stick, which were previously personalised and saved by the user through the **Reverse** button.  
Insert the memory stick in the appropriate USB port on the right side of the unit. Hold the **A/Rev** button down until a progress bar appears on the last line of the

## 5. Use

display. For each file, the unit will present torque, speed and reduction ratio values contained within the memory stick.

- ✓ **Reverse:** when the LED is on, the file will rotate clockwise if the pedal is held down or if the **Motor** button is pressed once. This function allows the file to be released once stopped as torque was reached within the root canal. **When the file rotates clockwise, the unit emits an intermittent acoustic signal.**

The torque, speed and reduction ratio personalised by the user for each file can be saved on the memory stick using the **Reverse** button.

Insert the memory stick in the appropriate USB port on the right side of the unit. Hold the **Reverse** button down until a progress bar appears on the last line of the display. For each file, the memory stick will contain the torque, speed and reduction ratio values **present on the unit at the time they are saved.**

**WARNING: all saving overwrites any previous saving.**


- ✓ **Red:** select the same reduction ratio as the reduction ratio of the contra-angle in use. If this does not occur, torque and rotation speed cannot be ensured by the unit. **WARNING:** if the reduction ratio is 1:1, the torque set is always equal to **Max** (no control).
- ✓ **Torque:** the torque scale is divided into 50 steps (from 0,2 Ncm to 5,1 Ncm to the file), plus a step relative to the lack of torque control by the unit (**Max torque**), in which the motor distributes all power available.

## 6. Use with EP-Fill (sold as an accessory)

GB

**ATR EndoPocket** is able to drive the **ATR EP-Fill** device. **ATR EP-Fill** is a handpiece that allows vertical obturation to be performed. Please consult the **ATR EP-Fill** instruction manual for any information concerning the device.

- ✓ Connect **ATR EP-Fill** and the pedal to **ATR EndoPocket**. Use the micromotor connector to connect **ATR EP-Fill**. Ensure the connection is well inserted.
- ✓ Once the **ATR EP-Fill** connector is connected to **ATR EndoPocket**, the following screen is displayed:

- ✓ Insert the plugger desired (see paragraph 4 Description of the ATR EP-Fill instruction manual).

The unit can be set to 3 different temperatures: 200°C, 250°C and 300°C, which can be selected using the Torque ▲▼ button.

**ATR EP-Fill** uses the pedal provided with **ATR EndoPocket** to activate the dental handpiece heating function. Pressing the pedal activates the heating function. Keeping the pedal pressed down allows **ATR EndoPocket** to reach the temperature set on the plugger and will keep it steady over time. During the period in which the temperature increases, the LED on the MOTOR button remains lit. Once the temperature set has been reached, **ATR EndoPocket** will emit a short acoustic signal. While the temperature is being kept steady over time, the LED on the MOTOR button will flash and be accompanied by a continuous acoustic signal.

**Keep the pedal pressed down to complete the procedure.**

If the pedal is released for a few seconds and then pressed down again, **ATR EndoPocket** will begin heating the plugger again starting at its current temperature; therefore the time needed to increase the temperature (MOTOR LED constantly lit) in this case can be quite brief or none at all.

GB

**6. Use with EP-Fill (sold as an accessory)**

For the entire time in which the pedal is pressed down and during which **ATR EndoPocket** manages **ATR EP-Fill** to reach the temperature set and keep it steady over time, the following will appear on the display:

**Warning:** do not hold the pedal down for more than 60 consecutive seconds. For 60 seconds of uninterrupted operation, wait 120 seconds before pressing the pedal again to activate further treatment.

**Allow 2 minutes for cooling before removing the plugger from the handpiece.**

**Warning:** the temperatures on the plugger can cause serious burns on human tissue. Only use ATR EP-Fill for filling while paying maximum attention to its use.

## 7. Cleaning, disinfection and sterilisation

GB

Clean the external surface of the handpiece and the relative cord using a cloth **slightly** dampened with water or **non-alcoholic disinfectant**.

- ✓ **The packaging is not airtight.**
- ✓ **Do not immerse the unit in ultrasound cleaners.**
- ✓ **Never, for any reason, spray liquid substances on the unit, especially on the monitor or in proximity of the electrical sockets.**

The micromotor, including the cable and connector, can be sterilised in an autoclave as long as the following instructions are followed:

- ✓ **Insert the motor safety cover.**
- ✓ **Place the micromotor, cable and connector in a bag suitable for sterilisation (refer to the recommendations by the autoclave manufacturer).**
- ✓ **Sterilise using a steam temperature of 121°C for 15 minutes.**
- ✓ **Do not use the micromotor again if its temperature is not below 40°C.**

Infection hazard if the following instructions are not complied with:

- ✓ **The sterilisation temperature must not exceed 135°C.**
- ✓ **Do not place the components in any kind of solution.**
- ✓ **Do not use ultrasound systems.**
- ✓ **Adequate maintenance and periodic operating and temperature checks of the autoclave are necessary.**

## 8. Maintenance

### *Servicing*

Check the micromotor cable at least once every 6 months. In the event the sheathing is worn, have it replaced by an authorised service centre.

**Do not lubricate the micromotor for any reason. After having lubricated the contra-angle, ensure the lubricant does not penetrate inside the micromotor.**

**Lubricant inside the micromotor strongly jeopardises unit operation and safety. In this condition, replacement of the micromotor cannot be considered an intervention covered by the warranty.**

### *Supplementary maintenance*

In the event supplementary maintenance is required, please contact **Advanced Technology Research S.r.l., Via S. Donato 1, 51100 - PISTOIA, ITALY ([www.atrdental.com](http://www.atrdental.com))** as manufacturer.

The manufacture is committed to providing the machine's electro-mechanical diagrams to personnel it deems qualified to perform repairs.

## 9. Specifications

GB

MANUFACTURER:	Advanced Technology Research S.r.l. PISTOIA – ITALY
MODEL:	<b>EndoPocket</b>
DIMENSIONS:	Open: 160x130x120 mm Closed: 160x130x60 mm
MATERIAL:	Packaging: bayblend Micromotor: aluminium
WEIGHT:	1.2 kg
VOLTAGE SUPPLY:	100-240 V
VOLTAGE FLUCTUATIONS:	MAX $\pm$ 10 %
FREQUENCY:	47-63 Hz
POWER RATING:	30 W
CURRENT RATING:	1.25 A
TORQUE VARIATION:	0,2-5,1 Ncm at file +/- 10%
SPEED VARIATION TO MOTOR SHAFT:	1600-12800 rpm +/- 10%
ELECTRIC SAFETY:	CLASS II
APPLIED COMPONENT TYPE:	BF (micromotor) 
LEVEL OF SAFETY IN PRESENCE OF INFLAMMABLE ANAESTHETIC MIXTURE OR OXYGEN:	NOT SUITABLE FOR USE IN PRESENCE OF INFLAMMABLE ANAESTHETIC MIXTURES OR OXYGEN
OPERATING MODE:	INTERMITTENT: 30 SEC. ON / 10 SEC. OFF
ENVIRONMENT CONDITIONS FOR USE:	Temperature +18°C to +40°C (+64°F to +104°F) Relative humidity: < 80%
PROTECTION AGAINST LIQUID PENETRATION:	ORDINARY
CLASSIFICATION AS MEDICAL DEVICE:	II, RULE IX 93/42/EEC
CENTRAL UNIT:	IP20
MOTOR:	IP65
PEDAL:	IP21
TRANSPORT AND STORING CONDITIONS:	+5°C to +65 C°(+41 to +149 °F); 20-95% RH NOT CONDENSING; 500-1060 hPa

Send your request to [product.service@atr dental.com](mailto:product.service@atr dental.com) to receive a list of authorised service centres. Customer Service direct line +39 0573 535846.

Upon receiving a written request, manufacturer Advanced Technology Research S.r.l., Via S. Donato 1, 51100-PISTOIA, ITALY will send a technical description, component list, calibration and testing instructions and all information necessary to allow qualified personnel to repair or maintain components it deems repairable.

## 10. Warranty Conditions

With this document, the manufacturer certifies the product has been properly assembled, top quality materials have been used, all necessary testing has been performed and product compliance with current legislation. The product is covered by a warranty period of 12 months, beginning the date of delivery to the final user, which must be attested to upon purchase by fiscal documentation issued by the retailer including the unit's serial number.

The consumer has the right to request warranty service during the period covered by the warranty only by sending a *written* notice no later than *2 months* from detecting the malfunction.

1. The warranty is limited to replacing or repairing components or pieces, which are the result of faulty production, with the exclusion of any expenses relating to labour, technicians' travel expenses, transportation and packaging, etc.
2. The warranty does not cover damage or malfunctions deriving from poor maintenance, incorrect power supply, negligence and incompetence. Likewise, the warranty does not cover malfunctions caused by the lack of servicing due to user negligence (see instruction manual).
3. This warranty does not include any direct or indirect compensation for damages of any nature to people or things as the result of the inefficiency of the equipment.
4. The warranty is automatically forfeited in the event the equipment is repaired, modified or tampered with by the purchaser or by unauthorised third parties.
5. For interventions covered by the warranty, the purchaser can only contact the retailer or the service centres indicated by the manufacturer or contact the manufacturer itself. Likewise the purchaser shall be responsible for ensuring the packaging is in perfect condition and suitable for safe transportation.
6. In the event of disputes concerning the application of the warranty or the quality or conditions of the equipment delivered, the purchaser cannot interrupt or delay the payment of the price or of the price instalments.
7. No compensation for damage can be requested by the purchaser due to inactive equipment.
8. The warranty explicitly excludes those defects:
  - caused by damage during transport;
  - that cannot be traced back to fabrication defects, but rather to normal wear and tear of the materials due to the proper use of the unit (specifically, calcification and wear and tear of the components subject to wear and tear)
  - caused by atmospheric events such as lightning and/or fire, and humidity and any other damage which cannot objectively be attributed to the manufacturer.

### RULES

- ✓ THE EFFECTIVENESS OF THE WARRANTY IS ENSURED BY THE PRESENCE OF THE FISCAL DOCUMENT ATTESTING THE PRODUCT PURCHASE DATE.
- ✓ THE WARRANTY DOES NOT COVER LABOUR, TRAVEL AND TRANSPORTATION EXPENSES, WHICH WILL ALWAYS BE AT THE EXPENSE OF THE PURCHASER.

## 11. Troubleshooting

GB

THE UNIT DOES NOT OPERATE PROPERLY	<p>Ensure the cable is properly plugged into the power socket and the network voltage corresponds to the voltage indicated on the power supply labelling.</p> <p>Press Function ▲ and Function ▼ together to restore factory parameters.</p>
THE MONITOR DOES NOT OPERATE PROPERLY	<p>If the monitor tends to become less bright, it means ATR EndoPocket is distributing the maximum power available. In this case, check the integrity of the contra-angle.</p>
THE MICROMOTOR DOES NOT START UP	<p>Ensure the motor connector is properly housed in the micromotor housing (red). Check the functional integrity of the contra-angle. Remove the contra-angle and set maximum speed, then start up the motor again. Calibrate without the contra-angle, then connect the contra-angle again and calibrate.</p>
THE PEDAL DOES NOT START UP THE MOTOR	<p>Try starting up the motor through the button on the keyboard. In the event the motor starts up, call the service centre to replace the pedal.</p>

**12. Declaration of Conformity**

Manufacturer: Advanced Technology Research A.T.R. S.r.l.  
Via S. Donato, 1 - 51100 Pistoia - Italy

Product: **EndoPocket**


Accessories included: Micromotor with protection cap for sterilisation, pedal, memory stick, external power supply and rubber micromotor support

The aforementioned products conform with the following European Directive:

Electromedical Equipment: Italian Leg. Decree 46/97 dated 24 February 1997 no. 46 and subsequent amendments, implementation of Directive "93/42/EEC".

April, 2007

ADVANCED TECHNOLOGY RESEARCH  
ATR srl

  
Jan Sicfert-PRESIDENT-

## 13. Accompanying Documents

GB

### *Electromagnetic emissions*

**ATR EndoPocket** is intended for use in the electromagnetic environment specified below. The customer or the user of the **ATR EndoPocket** should be used in such an environment.

RF emissions CISPR 11	Group 1	<b>ATR EndoPocket</b> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	<b>ATR EndoPocket</b> is suitable for use in all establishment including domestic establishments and those directly connected to the public low-voltage power supply.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuation/ flicker emissions IEC 61000-3-3	Conform	

## 13. Accompanying Documents

### *Electromagnetic immunity*

**ATR EndoPocket** is intended for use in the electromagnetic environment specified below. The **ATR EndoPocket** should only be used in such an environment.

Immunity test	Test level ICE 60601	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD)	$\pm 2,4,6\text{kV}$ contact $\pm 2,4,8\text{kV}$ air	$\pm 2,4,6\text{kV}$ $\pm 2,4,8\text{kV}$	Floor should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient / burst IEC 61000-4-4	$\pm 2\text{kV}$ for power supply lines $\pm 1\text{kV}$ for input / output lines	$\pm 2\text{kV}$ $\pm 1\text{kV}$	Mains power quality should be than of a typical commercial or hospital environment
Surge IEC 61000-4-5	$\pm 1\text{kV}$ differential mode $\pm 2\text{kV}$ cycle	$\pm 1\text{kV}$ $\pm 2\text{kV}$	Mains power quality should be than of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variation on power supply input lines.  IEC 61100-4-11	$<5\% U_t$ ( $>95\%$ dip in $U_t$ ) for 0,5 cycle  $40\% U_t$ ( $60\%$ dip in $U_t$ ) for 5 cycle  $70\% U_t$ ( $30\%$ dip in $U_t$ ) per 25 cycle  $<5\% U_t$ ( $>95\%$ dip in $U_t$ ) for 5 sec	$<5\% U_t$  $40\% U_t$  $70\% U_t$  $<5\% U_t$	Mains power quality should be than of a typical commercial or hospital environment. If the user of the <b>ATR EndoPocket</b> requires continued operation during power mains interruptions, it is recommended than the <b>ATR EndoPocket</b> be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic field IEC 61100-4-8	3A/m	N.A.	Power frequency magnetic fields should be at a typical commercial or hospital establishment level.

**13. Accompanying Documents**

**GB**

**ATR EndoPocket** is intended for use in the electromagnetic specified below. The customer or the user of the **ATR EndoPocket** should only be used in such an environment.

Immunity test	Test level IEC 60601	Compliance level	Electromagnetic environment - guidance
Conducted RF	3Vrms	3 V	<p>Portable and mobile RF communication equipment should be used no closer to any part of the <b>ATR EndoPocket</b>, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter</p> <p>Recommended separation distance</p>
IEC 61000-4-6	150 kHz to 80 MHz		$d = [3.5/V_1] \sqrt{P}$
Radiated RF	3V/m	3 V/m	$d = [3.5/E_1] \sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$
IEC 61000-4-3	80MHz to 2.5GHz		$d = [7/E_1] \sqrt{P} \quad 800\text{MHz to } 2.5\text{GHz}$ <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (w) according to the transmitter manufacturer and recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitter, as determined by an electromagnetic site survey <sup>a</sup> should be less than the compliance level in each frequency range <sup>b</sup>.</p> <p>Interference may occur in the proximity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>

### 13. Accompanying Documents

- a) Field strengths from fixed transmitters, such as base station for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted with accuracy. To assess the electromagnetic environment to fixed RF transmissions, an electromagnetic site survey should be considered if the measured field strength in the **ATR EndoPocket** is used exceeds the applicable RF compliance level above the **ATR EndoPocket** should be observed to verify normal operating. If abnormal performance is observed, additional measures may be necessary such as re-orienting or re-locating **ATR EndoPocket**.
- b) Over the frequency range 150kHz to 80MHz, field strengths should be less than [V] Vm.

#### **Recommended separation distance between portable and mobile RF communications equipment and the ATR EndoPocket**

The **ATR EndoPocket** is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the **ATR EndoPocket** can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **ATR EndoPocket** as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according in frequency of transmitter (m)		
	From 150kHz to 80MHz $d= 1,2\sqrt{P}$	From 80MHz to 800MHz $d= 1,2\sqrt{P}$	From 800MHz to 2,5GHz $d= 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

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